


HANAR STENT[®]


Gastro-Seal[™] Stent



Indicated for sealing leakage after sleeve gastrectomy

HANAR STENT®

Gastro-Seal™ Stent



“Endoscopic management post bariatric surgery leaks” presented in IFSO2016
 - Dr. Hany Shehab, Cairo University

Model Number	Stent		Delivery Device	
	Diameter Head-body-Head (mm)	Total Length (mm)	Diameter (mm)	Usable Length (mm)
ECBS-28-180-090	36-28-36	180	8	900
ECBS-28-210-090	36-28-36	210	8	900
ECBS-28-240-090	36-28-36	240	8	900